Catalog No. GIST01C

GIST-T1 Cell Line

(Human Cell Line)

May 7, 2019

Principle

Gastrointestinal stromal tumors (GISTs) are one of the submucosal tumors; occur in the stomach, the small intestine and the esophagus, unlike most gastrointestinal tumors. GISTs are considered to arise from the interstitial cells of Cajal, the pacemaker cells of the gut.

The GIST-T1 is a cell line derived from GISTs of the stomach in a Japanese woman and established by Takahiro Taguchi; associate professor, Graduate School of Integrated Arts and Sciences, Kochi-University, Kochi, Japan.

Warranty

Cosmo Bio warrants that product (cells) will be viable until the expiration date and is valid only if the product is stored and cultured according to the information indicated in this product data sheet. Cosmo Bio has optimized the cell culture media formulation which is ideal for the product. While other, unspecified cell culture media may also produce satisfactory results, a change in cell culture media or the absence of an additive(s) from the recommended cell culture media may affect recovery, growth and/or function of the product. If an alternative cell culture medium formulation is used to culture the product, the Cosmo Bio warranty for cell viability is no longer valid.

Precautions

- Because of cells derived from human tissue, please always wear gloves and safety glasses when working them.
- Remove the cryovial from the dry ice packaging and immediately place into liquid nitrogen storage until use.
- Based on the license agreement of Techno network Shikoku and Kochi University, GIST-T1 cell is prohibited to provide (distribution, lending, transfer, licensing, etc.) to a third party.

Components

Product Name	Size	Quantity	Storage Conditions	Expiration date
GIST-T1, cryopreserved	1.0x10 ⁶ cells/vial	1	Liquid Nitrogen vapor phase	6 months

*Shipping: dry ice

Optimized Culture Medium (selling separately)

Catalog No.	Product Name	Size	Quantity	Storage Conditions	Expiration date
GISTM	GIST-T1 Culture Medium	500 mL	1	-20°C Freezer	- Written on the bottle (stored at -20°C) 3months after thawing (stored at 4°C)
				4°C	
GISTMA	GIST-T1 Culture medium (without antibiotics)	500 mL	1	-20°C Freezer	
				4°C	

Culture Medium components: DMEM, FBS, etc.

General Information

Organism	<i>Homo sapiens</i> , human
Tissue	Stomach
Cultural Properties	Adherent
Biosafety Level	1
Gender	Female
Ethnicity	Asian
Virus Check	HIV-1(-), HTLV-1(-), HBV(-), HCV(-), T.pallidum(-)
Quality Check	Mycoplasma (-)

Materials required but not provided

- ✧ Variable volume pipettes
- ✧ Culture vessels
- ✧ 0.25% Trypsin
- ✧ HBSS or PBS(-)

Protocol

A) Thawing of Cells

- 1) Prepare a 100mm dish (**Note:** 100mm dish is recommended).
- 2) Warm culture medium to 37°C.
- 3) Prepare a conical tube (for 15mL) added 10mL of culture medium.
- 4) Carefully remove the cryovial from liquid nitrogen and thaw cells in a water bath at 37°C for 120 seconds.
- 5) Transfer the cryovial into a laminar flow hood. Before opening, wipe the outside of the vial with 70% ethanol.
- 6) Gently transfer the thawed cell suspension (1mL) into 10 mL of culture medium.
- 7) Transfer 1mL of culture medium in the same conical tube back to the cryovial and pour the contents back to 15mL conical tube.
- 8) Transfer the cell suspension to 100mm dish and incubate the cells in 37°C, 5% CO₂ incubator.
Note: Dilution and centrifugation of cells after thawing are not recommended since these actions are more harmful to the cells than the effect of residual DMSO in the culture.
- 9) Replace the medium with fresh pre-warmed culture medium every 2 to 3 days.

B) Subculturing Procedure

Note: Allow culture medium, HBSS(or PBS(-)), and 0.25% Trypsin to come to room temperature before use.

- 1) When the cells reach 70 -90% of confluent, they should be subcultured.
- 2) Aspirate the medium. Rinse the dish with 10mL of HBSS or PBS (-).
- 3) Add 1mL of 0.25% Trypsin, then incubate at 37°C for 4-6 minutes.
- 4) Add 10mL of culture medium and disperse the cells with gentle pipetting.
- 5) Transfer the cell suspension to conical tube and centrifuge at 200 ×g for 5 minutes at 4°C.
- 6) Aspirate the supernatant without disrupting the pellet and resuspend the cells in 10mL of culture medium.
- 7) Dilute the cell suspension by adding culture medium. A subcultivation ratio of 1:6 to 1:8 is recommended.
- 8) Transfer the cell suspension to new 100mm dish and Incubate the cells in 37°C, 5% CO₂ incubator.
- 9) Replace the medium with fresh pre-warmed culture medium every 2 to 3 days.
- 10) Culture the cells until the required density (70 -90% of confluent; Fig 1, C) is reached.

C) Freezing Procedure

- 1) Follow steps 1-6 [Subculturing Procedure] above.
- 2) Centrifuge the 50 mL tube at $\sim 200 \times g$ for 5 minutes at room temperature. Remove supernatant without disturbing the pelleted cells.
- 3) Resuspend cell pellet to $\sim 1 \times 10^6$ cells/mL with cold COS Banker Cell Freezing Medium (Cosmo Bio Co., Ltd. Cat No. KOJ-COS-CFM01) or other standard DMSO-based freeze media.
- 4) Dispense aliquots of this cell suspension into cryogenic storage vials. Keep vials in cell freezing container (e.g. CoolCell®, Biocision). Store vials at -80°C overnight.
- 5) Transfer the vials into gas phase above liquid nitrogen for long term storage.

Technical Information




Fig.1 Immunohistochemical and phase-contrast microscopic observation

A: Anti-CD34

B: Anti-c-kit

C: phase-contrast microscopic observation

COSMO BIO Co., LTD.

Inspiration for Life Science

References

- 1) Takahiro Taguchi, Hiroshi Sonobe, and Kazunari Yuri. et al. Conventional and Molecular Cytogenetic Characterization of a New Human Cell Line, GIST-T1, Established from Gastrointestinal Stromal Tumor. Lab Invest. 2002 May;82(5):663-5.
- 2) Watabe K, Sakamoto T, Kawazoe Y, Michikawa M, Miyamoto K, Yamamura T, Saya H, Araki N. Tissue culture methods to study neurological disorders: Establishment of immortalized Schwann cells from murine disease models. Neuropathology 2003;23:64-74.

COSMO BIO Co., LTD.

[JAPAN]

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: 81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

[Outside Japan]

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobioussa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioussa.com

For research use only. Not for clinical diagnosis.

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COSMO BIO CO., LTD.