


Assay for determining antimicrobial effects of solid polymer compounds on bacteria via incubation & CFU

Graphical Schematic of Procedure:



Materials:

- Sterile LB
- DPBS -CaCl₂, -MgCl₂ (Gibco Cat. No. 14190-094)
- 14 ml bacterial culture tubes with lids
- Falcon tubes 15 ml & 50 ml
- Pipette tips
- multichannel pipette
- 12-well plate with lid (round bottom)
- 96-well plate with lid (round bottom)
- LB agar plates
- Compounds for testing at desired concentration in LB
- Bacterial strain *E.coli*-AR3110 (or other) (-80°C stock)
- Eppendorf Bio Photometer
- Plastic cuvettes 1 ml
- Reservoir

Work under sterile conditions

Methods:

Day 1 → Streak out (single-colony method) bacterial strain (*E. coli*) on LB Agar plate as shown below (can use tip instead of toothpick) and allow to incubate overnight (ON) at 37°C.




Figure 1: Single colony streak out method

Day 2 → Prepare bacterial culture by inoculating 5mL LB with single colony from plate. Prepare negative control of only LB (3mL). Grow shaking 180-250 rpm at 37°C and allow to grow 6-16 hrs.

Confirm that the negative control LB is still clear and sterile. Measure OD₆₀₀ of bacterial culture and dilute to a starting culture of 0.01 using the $C_1V_1=C_2V_2$ formula. In a clean bench, add 500μL-750μL of 0.01 OD₆₀₀ diluted culture to the 12-well plate and place the sterile solid compound inside wells. Make sure to have proper controls: LB alone, LB + bacteria only, as well as additional positive or negative controls appropriate to the experiment. Cover plate with lid and place in incubator at 37°C with continuous shaking for 24 hours.

Day 3 → Look and take pictures of the results, use the plate directly for CFU assay.

Fill a 96-well plate with 180μL PBS (-,-) buffer using multichannel pipette via a reservoir. Add 20μL of incubated compound-bacterial solution that grew from previous day, to row A, each well has a different replicate/compound, and mix thoroughly. Change tips between compounds.

Serially dilute 20μL of mixture from row A to H, changing tips between rows and mixing thoroughly by pipetting up and down at least 10 times per well. May use multichannel pipette. Can throw away last 20μL from row H, to ensure same volume in each of the wells. (See diagram below).

Meanwhile dry LB-agar plates thoroughly on clean bench.

Mark LB-agar plates with a line using marker to be able to orient the lid.

Plate 20μL drops onto LB-agar starting from H to E, can use the same tip for each compound when increasing the dilution. Place lid on top of plate and trace where the droplets were placed onto the lid and label each drop. Carefully allow the drops to fully dry in clean bench, the colonies will become invisible once fully dry. See potential layout of 20μL drops below (Figure 3). Once dry cover and place upside down to grow ON in the incubator at 37°C.




Figure 2: 96-Well Plate Set-Up of Serial Dilution




Figure 3: Drop Colony Method for CFU. Each colony is 20 μ L. A standard LB plate can have between 16-25 20 μ L colonies.

Day 4 → Using colony counter and marker, count the colonies in each drop and write the numbers down and the corresponding dilution. Record the numbers in lab notebook. Analyze values in excel applying the following formula for the relationship between number of colonies, volume plated (20 μ L) and dilution factor:

$$\frac{CFU}{mL} = \frac{(\text{number of colonies} \times \text{dilution})}{\text{volume plated}}$$

Then graph the resulting CFU/mL as bar graphs using standard deviation values as error bars values.




Figure 4: LB Agar plate after colonies have formed and counted with marker