

## Standard Operation Procedure (SOP)

### Inactivating Genetically Modified Organisms (GMOs) with Formaldehyde Solution

The SOP refers to an inactivation of the sample material described below and under the conditions of use here:

#### Preconditions and precautions:

- Get properly introduced and trained to work with bacteria or viruses.
- Wear two pairs of gloves nitrile or latex gloves over each other, sleeve protectors and a clean lab coat. If gloves get possibly contaminated during the procedure, change them immediately.
- Perform all step of the procedure inside a sterile work bench (Class II) until samples are covered or resuspended in the formaldehyde solution. Preparation and removal of the formaldehyde solution have to take place in a chemical fume hood.
- Disinfect sterile work bench (Class II) with Meliseptol before and after finishing work.
- The formaldehyde solution need to be freshly prepared; in case it is older than 3 months, prepare a fresh solution according to the protocol below and document in the formaldehyde preparation log book and on the respective formaldehyde bottle!

#### Reagents:

- Meliseptol
- 37 % formaldehyde solution: ROTIPURAN 37 %, p.a., ACS, article no. 4979.1, Roth
- 1x PBS<sup>-</sup>, pH 7.4 without calcium and magnesium (1000 mL): Dilute 100 ml 10x PBS (ROTI®Fair 10x PBS 7.4, article no. 1105.1, Roth in 1000 mL MilliQ water) in 900 ml of MilliQ water.
- 1x DPBS<sup>++</sup>, pH 7.4 with Calcium and magnesium (Gibco)

#### Equipment and materials:

- Tubes or plates with lid
- Parafilm
- MilliQ water
- Vacuboy acceleration system & sterile, one-time usable plastic Pasteur pipettes
- Pipette boy & serological pipettes: 5 ml, 10m ml and 1000 ml
- Micropipette & micropipette tips: 10 µl, 100 µl, 10000 µl
- Centrifuge with aerosol tight rotor and aerosol tight rotor lids
- Sterile work bench class II
- Chemical fume hood
- 37 °C incubator

#### Virus or bacteria sample material:

1. Recombinant human herpesvirus 1 (**HSV-1**) of risk group 2 (produced in S2 work 92/14-2)
2. Recombinant influenza A virus (**IAV**) of risk group 2 (produced in S2 work 92/14-3)
3. Genetically modified **E. coli** Bacteria of risk group 1 such as K12-derivatives (produced or used e.g. in S1 work n° 10 of the genetic engineering annex 92/14)
4. as well as **eukaryotic cell cultures** of risk group 1, which may have been infected with the above-mentioned organisms.

#### Application conditions:

1. Incubation of the sample material such as **HSV-1**, **IAV** and **E. coli K12 derivatives** (see above) in a **4 % formaldehyde solution** for **30 minutes** at **room temperature** (20-25°C).
2. Incubation of cell-free recombinant influenza A viruses (IAV) in a minimum **0.015 % formaldehyde solution** for **16 hours** at **37°C**.

### Preparation of formaldehyde solutions:

- For 4 % formaldehyde solution (500 mL):  
Add 54 mL 37 % formaldehyde solution to 446 mL of 1x PBS<sup>-</sup>.
- For 0.015 % formaldehyde solution (500 mL):  
Add 2 mL freshly prepared 4 % formaldehyde solution (see above) to 498 mL of 1x PBS<sup>-</sup>.

*Document the preparation of formaldehyde solutions in the log book in lab 114.*

### Procedure for formaldehyde inactivation:

- Transfer closed tube or well plate containing virus or bacteria to a sterile work bench Class II
- Wash the virus or bacteria sample with 1x DPBS<sup>++</sup>, pH7.4 and add 4 % Formaldehyde solution (not older than 3 month); in detail:
  - o *For virus or bacteria grown on/in cells attached to plates,*
    - discard supernatant on cells with a sterile Pasteur pipette and the Vacuboy (take care not to scrape cells off of the well surface);
    - wash cells with 1x DPBS<sup>++</sup>, pH7.4;
    - Discard supernatant on cells with a sterile Pasteur pipette and the Vacuboy (take care not to scrape cells off of the well surface);
    - cover cells completely with 4 % Formaldehyde solution.
  - o *For virus or bacteria in solution,*
    - centrifuge bacteria or viruses in solution in tubes/plates with closed lids and using an aerosol tight rotor with aerosol tight rotor lid applied;
    - resuspend virus or bacteria pellet in DPBS<sup>++</sup>, pH 7.4;
    - centrifuge bacteria or viruses in solution in tubes/plates with closed lids and using an aerosol tight rotor with aerosol tight rotor lid applied;
    - resuspend virus or bacteria pellet in 4 % Formaldehyde solution or only for cell-free IAV samples only in minimum 0.015 % formaldehyde solution.
- Close Formaldehyde container immediately after use as well as the tube or plate containing the virus or bacteria sample in formaldehyde solution.
- Clean well closed plates or tubes containing the virus or bacteria sample with Meliseptol, take out of the sterile work bench and incubate for proper inactivation with tightly closed lid:
  - o *For samples in 4 % formaldehyde solution,*  
incubate for 30 min (do not incubate less, you can incubate for longer) at room temperature (20-25°C) in the chemical fume hood or
  - o *For cell-free IAV samples only in minimum 0.015 % formaldehyde solution,*  
incubate for 16 hours at 37°C in a well closed and with parafilm sealed container in at 37 °C in an incubator.
- After the proper inactivation time, transfer the samples with tightly closed lids to a chemical fume hood and discard formaldehyde solution into the formaldehyde waste container:
  - o *For virus or bacteria grown on/in cells attached to plates by accelerating the formaldehyde with a pipette boy and a serological pipette followed by adding MilliQ water or 1x PBS<sup>-</sup>, pH 7.4.*
  - o *For virus or bacteria in solution, by first centrifuging the tightly closed tubes/plates using an aerosol tight rotor with aerosol tight rotor lid applied and then accelerating the formaldehyde with a pipette boy and a serological pipette;*  
resuspend pellet in MilliQ water or 1x PBS<sup>-</sup>, pH 7.4.
- Store properly labelled sample at 4 °C for further applications.