Area of application: Freie Universität Berlin Altensteinstraße 23a



Operating Instructions – Heating stirrer According to Occupational Health and Safety Regulation		
Work	space: SupraFAB Work s	station: Room 114/115.7./114.3
Hazards for people and the environment		
EX	Danger of skin burns Danger due to splashing of the reaction mixture by the magne	tic fish.
Protective measures and rules of behaviour		
	The system must be tested for operational stability before use Open containers must be protected against accidental splashir If necessary, wear suitable protective gloves and eyeglasses.	
Behaviour in case of danger		
1	Inform the person responsible for the area immediately in the machine! In case of fire: Keep calm! Call the fire department on 112! Inform the central control re Warn employees! Inform the person responsible for the area! Fight the fire with a powder, CO2 or foam extinguisher! Pay at	oom on <mark>55112</mark> !
First aid		
Rescue service: 112 Rinse burnt main areas, allow to bleed out if possible and spray with disinfectant immediately, add disinfectant if necessary and leave to work for at least 30 minutes. If necessary, alert first aiders, emergency services or a doctor. Inform supervisor. Every injury, no matter how minor, must be entered in the first aid book		
Service		
	Maintenance and repair work may only be carried out with the written permission of the laboratory manager. The function of the heating stirrer must be checked before commissioning. Only have maintenance and repairs carried out by a specialist.	
	Only spare parts that correspond to the original parts in terms of material and design may be used for maintenance.	
Observe the manufacturer's maintenance schedule.		
For further questions, please contact the Occupational Safety department Tel.: -54495 Date and Signature of the division manager:		